# VarioskanLUX

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
The screenshot shows a software interface with a top navigation bar and a main workspace. The top bar includes a '結果' (Results) tab and a '計算の追加' (Add Calculation) button. Below this, there are several icons for different calculation types: 'ブランク減算' (Blank subtraction), '平均、SD、CV%' (Average, SD, CV%), '基本計算' (Basic calculation), '1:X 希釈率' (1:X Dilution rate), '% 正規化' (% Normalization), '光路長補正' (Path length correction), '標準曲線' (Standard curve), '用量反応' (Dose response), '平行線分析' (Parallel line analysis), '酵素のカイネティクス' (Enzyme kinetics), and 'カイネ' (Kinetics). The 'ブランク減算' and '平均、SD、CV%' icons are highlighted with a red box. The main workspace shows a '新しいセッション\*' (New session\*) window with a tree view on the left and a data table on the right. The tree view has '結果' (Results) expanded, showing '吸光度 1' (Absorbance 1) and 'ブランク減算 1' (Blank subtraction 1). The '平均、SD、CV%' icon is also highlighted with a red box. The data table has columns for '1', '2', and '3' and rows for 'A', 'B', and 'C'. The table content is as follows:

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The screenshot shows the SkanIt Software 7.0.2 RE for Microplate Reader interface. The title bar reads "新しいセッション\* - SkanIt Software 7.0.2 RE for Microplate R". The main window is titled "新しいセッション\*" and contains a left sidebar with a tree view and a main workspace.

**Protocol Setup:**

- The "プロトコル" (Protocol) step is selected in the tree view.
- Under "プロトコル", two sub-steps are added: "吸光度 1" (Absorbance 1) and "吸光度 2" (Absorbance 2).
- In the main workspace, the "波長" (Wavelength) field is set to "620" nm.
- The "マルチ波長" (Multi-wavelength) checkbox is unchecked.
- The "光路長補正" (Optical path length correction) checkbox is unchecked.
- The "高度なパラメータ" (Advanced parameters) section is expanded.
- The "開始" (Start) button is highlighted with a red box.

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